



OccuSwitch

LRH1070/00 SENSR SURFACE BOX

Ceiling mounting box for the OccuSwitch LRM1070/1080


Ceiling mounting box for the OccuSwitch LRM1070/1080

Product data

Product Data	
Full product code	871155973143899
Order product name	LRH1070/00 SENSR SURFACE BOX
EAN/UPC - Product	8711559731438
Order code	913700328003
Numerator - Quantity Per Pack	1

Numerator - Packs per outer box	42
Material Nr. (12NC)	913700328003
Net Weight (Piece)	0.100 kg

Dimensional drawing



LRH1070/00 SENSR SURFACE BOX