



Group Cabinet Controls

LCU7591 Amplight Leak Coil

Leak Coil / current-sensor to be used in combination with the Current to measure the leakage current


Leak Coil / current-sensor to be used in combination with the Current to measure the leakage current

Product data

Product data		Product data	
Full product code	872790094774800	Numerator - Packs per outer box	18
Order product name	LCU7591 Amplight Leak Coil	Material Nr. (12NC)	913700341903
EAN/UPC - Product	8727900947748	Net Weight (Piece)	0.425 kg
Order code	913700341903		
Numerator - Quantity Per Pack	1		

Group Cabinet Controls

Dimensional drawing



LCU7591 Amplight Leak Coil

