



Fortimo LED strip



Fortimo LED Strip 2ft 2200lm 930 HV4 +

Fortimo LED Strip CED BioUp is a carrier for Signify's BioUp dynamic technology. This technology provides a dynamic CCT with high MDER functionality. CCT tunable from 2700K to 5000K. The modules are designed to be paired with Xitanium FlexTune drivers to enable an easy set-up and programming.

Product data

General Information	
Lighting Technology	LED
Light Technical	
Luminous Flux	2,073 lm
Life to 70% Lumen Maintenance	70,000 hour(s)
Correlated Color Temperature (Nom)	3000 K
Luminous Efficacy (rated) (Nom)	168 lm/W
Color Consistency	<3
Color rendering index (CRI)	90
Temperature	
T-Case Lifetime (Max)	70 °C
Mechanical and Housing	
Net Weight (Piece)	0.040 kg
Approval and Application	
Energy Efficiency Class	D

EPREL Registration Number	766628
Remarks	Please download the Lighting - Product Data Sheet for more information and ordering options
Product Data	
Order product name	Fortimo LED Strip 2ft 2200lm 930 HV4 +
Full product name	Fortimo LED Strip 2ft 2200lm 930 HV4 +
Full product code	871869969426500
Order code	929002111406
Material Nr. (12NC)	929002111406
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718699694265
Numerator - Packs per outer box	168
EAN/UPC - Case	8718699694272