



# Fortimo LED strip



## Fortimo LED Strip 0.5ft 490lm 940 FC HV5

Fortimo LED Strip CED BioUp is a carrier for Signify's BioUp dynamic technology. This technology provides a dynamic CCT with high MDER functionality. CCT tunable from 2700K to 5000K. The modules are designed to be paired with Xitanium FlexTune drivers to enable an easy set-up and programming.

### Product data

General Information	
Lighting Technology	LED
Light Technical	
Luminous Flux	467 lm
Life to 70% Lumen Maintenance	60,000 hour(s)
Color Consistency Initial	3
Correlated Color Temperature (Nom)	4000 K
Luminous Efficacy (rated) (Nom)	164 lm/W
Color Consistency	<3
Color rendering index (CRI)	90
Temperature	
T-Case Lifetime (Max)	70 °C
T-Case Maximum (Max)	80 °C
Mechanical and Housing	
Overall length	140 mm
Overall width	20 mm
Overall height	6.1 mm
Dimensions (Height x Width x Depth)	6 x 20 x 140 mm

Net Weight (Piece)	0.011 kg
Approval and Application	
Energy Efficiency Class	C
EPREL Registration Number	778292
Remarks	Please download the Lighting - Product Data Sheet for more information and ordering options
Product Data	
Order product name	Fortimo LED Strip 0.5ft 490lm 940 FC HV5
Full product name	Fortimo LED Strip 0.5ft 490lm 940 FC HV5
Full product code	871869972148000
Order code	929002119406
Material Nr. (12NC)	929002119406
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718699721480
Numerator - Packs per outer box	56
EAN/UPC - Case	8718699721497