



Fortimo LED strip



Fortimo LED Strip 1400mm 5500lm 930 LV5F

Fortimo LED Strip CED BioUp is a carrier for Signify's BioUp dynamic technology. This technology provides a dynamic CCT with high MDER functionality. CCT tunable from 2700K to 5000K. The modules are designed to be paired with Xitanium FlexTune drivers to enable an easy set-up and programming.

Product data

General Information	
Lighting Technology	LED
Light Technical	
Luminous Flux	4,925 lm
Life to 70% Lumen Maintenance	60,000 hour(s)
Color Consistency Initial	3
Correlated Color Temperature (Nom)	3000 K
Luminous Efficacy (rated) (Nom)	152 lm/W
Color Consistency	<3
Color rendering index (CRI)	90
Temperature	
T-Case Lifetime (Max)	70 °C
T-Case Maximum (Max)	80 °C
Mechanical and Housing	
Overall length	1,399.2 mm
Overall width	20 mm
Overall height	6.1 mm
Dimensions (Height x Width x Depth)	6 x 20 x 1399 mm

Net Weight (Piece)	0.080 kg
Approval and Application	
Energy Efficiency Class	D
EPREL Registration Number	778367
Remarks	Please download the Lighting - Product Data Sheet for more information and ordering options
Product Data	
Order product name	Fortimo LED Strip 1400mm 5500lm 930 LV5F
Full product name	Fortimo LED Strip 1400mm 5500lm 930 LV5F
Full product code	871869975959900
Order code	929002177806
Material Nr. (12NC)	929002177806
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718699759599
Numerator - Packs per outer box	96
EAN/UPC - Case	8718699759605