



Fortimo LED strip



Fortimo LED Strip 0.25ft245lm 930 HV5F

Fortimo LED Strip CED BioUp is a carrier for Signify's BioUp dynamic technology. This technology provides a dynamic CCT with high MDER functionality. CCT tunable from 2700K to 5000K. The modules are designed to be paired with Xitanium FlexTune drivers to enable an easy set-up and programming.

Product data

General Information	
Lighting Technology	LED
Light Technical	
Luminous Flux	222 lm
Life to 70% Lumen Maintenance	70,000 hour(s)
Color Consistency Initial	3
Correlated Color Temperature (Nom)	3000 K
Luminous Efficacy (rated) (Nom)	155 lm/W
Color rendering index (CRI)	90
Temperature	
T-Case Lifetime (Max)	70 °C
T-Case Maximum (Max)	80 °C
Mechanical and Housing	
Overall length	70 mm
Overall width	20 mm
Overall height	6.1 mm
Dimensions (Height x Width x Depth)	6 x 20 x 70 mm

Net Weight (Piece)	0.010 kg
Approval and Application	
Energy Efficiency Class	D
Approval Marks	CB Certificate CE marking ENEC certificate
EPREL Registration Number	1227093
Product Data	
Order product name	Fortimo LED Strip 0.25ft245lm 930 HV5F
Full product name	Fortimo LED Strip 0.25ft245lm 930 HV5F
Full product code	871951433190700
Order code	929002919006
Material Nr. (12NC)	929002919006
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8719514331907
Numerator - Packs per outer box	56
EAN/UPC - Case	8719514331914