



Fortimo LED strip



Fortimo LED Strip 0.5ft MF 940 FC LV6

Fortimo LED Strip CED BioUp is a carrier for Signify's BioUp dynamic technology. This technology provides a dynamic CCT with high MDER functionality. CCT tunable from 2700K to 5000K. The modules are designed to be paired with Xitanium FlexTune drivers to enable an easy set-up and programming.

Product data

General Information		Dimensions (Height x Width x Depth)	
Lighting Technology	LED		4 x 20 x 140 mm
		Net Weight (Piece)	0.009 kg
Light Technical		Approval and Application	
Luminous Flux	550 lm	Energy Efficiency Class	D
Color Consistency Initial	3	Approval Marks	CE marking ENEC certificate CB Certificate
Correlated Color Temperature (Nom)	4000 K	EPREL Registration Number	1413389
Luminous Efficacy (rated) (Nom)	173 lm/W		
Color rendering index (CRI)	90	Product Data	
Operating and Electrical		Order product name	Fortimo LED Strip 0.5ft MF 940 FC LV6
Power Consumption	3.173 W	Full product name	Fortimo LED Strip 0.5ft MF 940 FC LV6
Temperature		Full product code	871951448666900
T-Case Lifetime (Max)	75 °C	Order code	929002940306
T-Case Maximum (Max)	85 °C	Material Nr. (12NC)	929002940306
Mechanical and Housing		Numerator - Quantity Per Pack	1
Overall length	140 mm	EAN/UPC - Product/Case	8719514486669
Overall width	20 mm	Numerator - Packs per outer box	66
Overall height	4.3 mm	EAN/UPC - Case	8719514486676