



Fortimo LED strip



Fortimo LED Strip 2ft MF 840 BC HV6

Fortimo LED Strip CED BioUp is a carrier for Signify's BioUp dynamic technology. This technology provides a dynamic CCT with high MDER functionality. CCT tunable from 2700K to 5000K. The modules are designed to be paired with Xitanium FlexTune drivers to enable an easy set-up and programming.

Product data

General Information	
Lighting Technology	LED
Light Technical	
Luminous Flux	2,200 lm
Color Consistency Initial	3
Correlated Color Temperature (Nom)	4000 K
Luminous Efficacy (rated) (Nom)	186 lm/W
Color rendering index (CRI)	80
Operating and Electrical	
Power Consumption	11.8014 W
Temperature	
T-Case Lifetime (Max)	75 °C
T-Case Maximum (Max)	85 °C
Mechanical and Housing	
Overall length	559.2 mm
Overall width	20 mm
Overall height	7.4 mm

Dimensions (Height x Width x Depth)	7 x 20 x 559 mm
Net Weight (Piece)	0.035 kg
Approval and Application	
Energy Efficiency Class	C
Approval Marks	CE marking ENEC certificate CB Certificate
EPREL Registration Number	1413442
Product Data	
Order product name	Fortimo LED Strip 2ft MF 840 BC HV6
Full product name	Fortimo LED Strip 2ft MF 840 BC HV6
Full product code	871951448912700
Order code	929002945606
Material Nr. (12NC)	929002945606
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8719514489127
Numerator - Packs per outer box	110
EAN/UPC - Case	8719514489134