



## LLC7871/00 CT NODE 20MM3 DC EU4VF DG

### Product data

Operating and Electrical		Order product name	LLC7871/00 CT NODE 20MM3 DC EU4VF DG
Input Current (Nom)	0.025 A	Order code	913701043603
Temperature		Numerator - Quantity Per Pack	1
Ambient temperature range	-40 to +70 °C	Numerator - Packs per outer box	45
Approval and Application		Material number (12NC)	913701043603
CE mark	CE mark	Full product name	LLC7871/00 CT NODE 20MM3 DC EU4VF DG
EU RoHS compliant	Yes	EAN/UPC - Case	8718699708924
Product Data			
Full product code	871869970891700		