



## LLC7856/00 CT NODE ZHAGA DC MS EU4VF LG

### Product data

Operating and Electrical		Order product name	LLC7856/00 CT NODE ZHAGA DC MS EU4VF LG
Input Current (Nom)	0.025 A	Order code	913701059403
Temperature		Numerator - Quantity Per Pack	1
Ambient temperature range	-40 to +70 °C	Numerator - Packs per outer box	45
Approval and Application		Material number (12NC)	913701059403
CE mark	CE mark	Full product name	LLC7856/00 CT NODE ZHAGA DC MS EU4VF LG
EU RoHS compliant	Yes	EAN/UPC - Case	8719514417304
Product Data			
Full product code	871951441729800		