



# Fortimo LED strip



## Fortimo LED Strip 1ft 1100lm 930 HV4 +

Fortimo LED Strip CED BioUp is a carrier for Signify's BioUp dynamic technology. This technology provides a dynamic CCT with high MDER functionality. CCT tunable from 2700K to 5000K. The modules are designed to be paired with Xitanium FlexTune drivers to enable an easy set-up and programming.

### Product data

General Information	
Remarks	Please download the Lighting - Product Data Sheet for more information and ordering options
Lighting Technology	LED
Light Technical	
Luminous Flux	1,036 lumen
Life to 70% Lumen Maintenance	70,000 hour(s)
Correlated Color Temperature (Nom)	3000 K
Luminous Efficacy (rated) (Nom)	168 lm/W
Color Consistency	<3
Color rendering index (CRI)	90
Temperature	
T-Case Lifetime (Max)	70 °C

Approval and Application	
Energy Efficiency Class	D
EPREL Registration Number	766625
Product Data	
Full product code	871869969418000
Order product name	Fortimo LED Strip 1ft 1100lm 930 HV4 +
Order code	929002111006
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	168
Material Nr. (12NC)	929002111006
Full product name	Fortimo LED Strip 1ft 1100lm 930 HV4 +
EAN/UPC - Case	8718699694197