



# Fortimo LED strip



## Fortimo LED Strip 2ft 2200lm 830 FC LV5

Fortimo LED Strip CED BioUp is a carrier for Signify's BioUp dynamic technology. This technology provides a dynamic CCT with high MDER functionality. CCT tunable from 2700K to 5000K. The modules are designed to be paired with Xitanium FlexTune drivers to enable an easy set-up and programming.

### Product data

General Information	
Remarks	Please download the Lighting - Product Data Sheet for more information and ordering options
Lighting Technology	LED
Light Technical	
Luminous Flux	2,000 lm
Life to 70% Lumen Maintenance	70,000 hour(s)
Color Consistency Initial	3
Correlated Color Temperature (Nom)	3000 K
Luminous Efficacy (rated) (Nom)	165 lm/W
Color Consistency	<3
Color rendering index (CRI)	80
Temperature	
T-Case Lifetime (Max)	70 °C
T-Case Maximum (Max)	80 °C
Mechanical and Housing	
Overall length	559.2 mm

Overall width	20 mm
Overall height	7.35 mm
Dimensions (Height x Width x Depth)	7 x 20 x 559 mm

Approval and Application	
Energy Efficiency Class	C
EPREL Registration Number	778332

Product Data	
Order product name	Fortimo LED Strip 2ft 2200lm 830 FC LV5
Full product name	Fortimo LED Strip 2ft 2200lm 830 FC LV5
Full product code	871869971729200
Order code	929002125206
Material Nr. (12NC)	929002125206
Numerator - Quantity Per Pack	1
EAN/UPC - Product/Case	8718699717292
Numerator - Packs per outer box	84
EAN/UPC - Case	8718699717308