



# Fortimo LED strip



## Fortimo LED Strip 1400mm 5500lm 840 HV5F

Fortimo LED Strip CED BioUp is a carrier for Signify's BioUp dynamic technology. This technology provides a dynamic CCT with high MDER functionality. CCT tunable from 2700K to 5000K. The modules are designed to be paired with Xitanium FlexTune drivers to enable an easy set-up and programming.

### Product data

General Information	
Remarks	Please download the Lighting - Product Data Sheet for more information and ordering options
Lighting Technology	LED
Light Technical	
Luminous Flux	5,315 lumen
Life to 70% Lumen Maintenance	70,000 hour(s)
Color Consistency Initial	3
Correlated Color Temperature (Nom)	4000 K
Luminous Efficacy (rated) (Nom)	174 lm/W
Color Consistency	<3
Color rendering index (CRI)	80
Temperature	
T-Case Lifetime (Max)	70 °C
T-Case Maximum (Max)	80 °C
Mechanical and Housing	
Overall length	1,399.2 mm

Overall width	20 mm
Overall height	6.1 mm
Dimensions (Height x Width x Depth)	6 x 20 x 1399 mm
Approval and Application	
Energy Efficiency Class	C
EPREL Registration Number	778358
Product Data	
Full product code	871869976078600
Order product name	Fortimo LED Strip 1400mm 5500lm 840 HV5F
Order code	929002176906
Numerator - Quantity Per Pack	1
Numerator - Packs per outer box	96
Material Nr. (12NC)	929002176906
Full product name	Fortimo LED Strip 1400mm 5500lm 840 HV5F
EAN/UPC - Case	8718699760793
EAN/UPC - Product/Case	8718699760786