

PHILIPS

Eco passport

LED Linear MOD 20W 865 2PK

MyCare LED Circular

MyCare LED Ceiling Module integrates a LED light source into a module. It is specially designed to replace fluorescent circular tubes working on a G10q fitting. The module offers significant energy savings and a long lifetime.



Our eco passport

Signify leads the way to brighter lives and a better world. As part of the design process, we innovate to reduce our environmental impact by taking climate action and by focusing on circularity and bring wider benefits to the world we operate in via social aspects as food availability, health & wellbeing, and safety & security.

We know this is important to you too. So we have developed an eco passport explaining the environmental and social performance of our products in eight Sustainable Focal Areas, developed through our Sustainable Design process since 1994.



Energy

- **Dimmable:** No
- **Luminous Efficacy (rated) (Nom):** 115 lm/W
- **Power (Nom):** 20 W



Circularity

- **Nominal Lifetime (Nom):** 20000 h



Packaging

- **SAP Weight Paper Unit (Case):** 2.500 kg



Substances

- **EU RoHS compliant:** Yes



Weight & Materials

- **Net Weight (Piece):** 0.125 kg

Further reading at:
<https://www.signify.com/global/sustainability/sustainable-lighting>

© 2023 Signify Holding. All rights reserved.