



# Cabinet Controls

## LCU7591 Amplight Leak Coil

Leak Coil / current-sensor to be used in combination with the Current to measure the leakage current


Leak Coil / current-sensor to be used in combination with the Current to measure the leakage current

### Product data

Product Data	
Full product code	872790094774800
Order product name	LCU7591 Amplight Leak Coil
EAN/UPC - Product	8727900947748
Order code	913700341903
Numerator - Quantity Per Pack	1

Numerator - Packs per outer box	18
Material Nr. (12NC)	913700341903
Net Weight (Piece)	0.425 kg

### Dimensional drawing



LCU7591 Amplight Leak Coil