



# MASTER LEDtube EM/ Mains T8

## MASTER LEDtube Starter EMP GenII OP

The Philips MASTER LEDtube integrates a LED light source into a traditional fluorescent form factor. Its unique design creates a perfectly uniform visual appearance which cannot be distinguished from traditional fluorescent. For those that are looking for value for money within limited budget and re-lamping efforts for better light effect and lifetime.

### Warnings and Safety

- NOTES: The overall energy efficiency and light distribution of any installation that uses these lamps are determined by the design of the installation.

### Product data

General Information		Order code	
EU RoHS compliant	Yes	Order code	929001801432
Product Data		Numerator - Quantity Per Pack	1
Full product code	871869681965400	Numerator - Packs per outer box	40
Order product name	MASTER LEDtube Starter EMP GenII OP	Material Nr. (12NC)	929001801432
EAN/UPC - Product	8718696819654	Net Weight (Piece)	0.004 kg

# MASTER LEDtube EM/Mains T8

## Dimensional drawing



Product	D1	D2	C1
MASTER LEDtube Starter EMP GenII OP	3 mm	21.5 mm	34.5 mm

Dummy starter for cost down version

